

Deluxe Sentra Full Reclining Wheelchair

HCPCS: K0001 OR K0002+E1226+E0990 OR K0195,K0001 OR K0002+E1226+E0990 OR K0195



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Durable, heavy-gauge naugahyde upholstery
- Carbon Steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Composite Mag-style wheels with chrome hand rims are lightweight and maintenance free
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels set back on frame prevent tipping
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Rear anti-tippers standard
- Comes standard with swing-away elevating legrests
- Aluminum foot plates are attractive, crack-proof and lightweight
- Front caster forks are adjustable in two positions
- Carry pocket standard
- Headrest extension with cushioned head immobilizer standard
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD22RBDDA	Sentra,22"Recl.DetDskArm,1/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	450 lbs
Seat To Floor Height	21"
Carton Shipping Weight	72 lbs
Back Height	23"J33"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
STDS1074	6"Hemi Caster,SentraHD,1ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS802	Anti-Tipprs,w/Wheels,Adj,1pr/bx, BX
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS

ACCESSORIES	Description
STDS831	Heel Loops, 1pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS834	OverheadAnti- Theft Bar,Single,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS6008-1	Carry Bag,W/C,Oxygen,1/ea, EA
PARTS	Description
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, EA
STDS4035R	PushToLockBrakeSentraRecliner,Right,1/ea, EA
STDS4035L	PushToLockBrakeSentraRecliner,Left,1/ea, EA
STDS1001	LTWT 24" wheel w/alum Handrim, 1/ea, EA
STDS818	Anti-Tippers,rear,SentraRecliner,1pr/bx, BX
STDS004010	Handrim,24"Wheel,Chrome/1ea, EA
STDS3003	Caster Fork,Whch,1/ea, EA
STDS004011	Handrim,Alloy24"Wheel/1ea., EA
STDS2324BLK	Seat Uphol 16"x16"Blk,1/bx, BX
STDS2A2942	Seat Uphol.22"SentraRecl.1/ea, EA
STDSFLVB	UniversalWheelChairVinylFullArmPad,1/ea, EA
STDS4014R	Piston,Sentra Recl.Rt.1/ea, EA
STDS4014L	Piston,Sentra Recl.Left,1/ea, EA
STDS1000	Wheel,24"Chr.Hand Rl.Sentra1/ea, EA
STDS1006	24" Wheel,Black Handrim,1/ea, EA
STDS1005-1	RearAxle Bolt asy, ViperPlus Recl, 1/pkg, PKG
STDS1012HD	Wheel Bearing,24"HD W/C,1/ea, EA
STDSDDASL	DeskArmw/Pads,Sentra,Lft.1/ea, EA
STDSDDASR	Desk Arm w/Pads,Sentra,Rt.1/ea, EA
STDSADDASPL	Desk Arm,Adj.Ht.Left,1/ea, EA