

Lightweight Expedition Aluminum Transport Chair

HCPCS: E1038



FIGURE A




FIGURE B

- Weighs only 19 lbs
- Available in blue or red frames
- Folds flat for easy transport (Figure A)
- Handy lever on chair back makes it easy and convenient to fold back down (Figure B)
- Comes with seat belt
- Companion activated wheel locks
- 12" rear "flat-free" wheels

PRODUCTS	Description
EXP19LTBL	Expedition X-Light Blue TransChair,1/ea
EXP19LTRD	Expedition X-Light Red TransChair,1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	19" (W) x 16" (D)
Weight Capacity	300 lbs
Seat To Floor Height	19"
Carton Shipping Weight	24 lbs
Overall Width	22" (open) 10.75" (closed)
Weight	19 lbs
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

PARTS	Description
STDS1070	Seat Rail Guide, Removable Arm, 1/ea, EA
STDS871	Hand Grip,Sentra,1/ea, EA
STDS4020	Tipping Sleeve 1/ea, EA
STDS002001	Arm Pad Screw,1/ea, EA
STDSFLVB	UniversalWheelChairVinylFullArmPad,1/ea, EA
STDS5D5505	Wheel Brake,Left,EXP19LT,1ea, EA
STDS5D5506	Wheel Brake,Right,EXP19LT,1ea, EA
STDS5D5507	Rear Axle Bearings,EXP19LT,1pair, PR
STDS5D5508	Rear Hub Cap,EXP19LT,1ea, EA
STDS5D5509	Rear Axle Bolt with Sleeve,EXP19LT,1ea, EA

PARTS	Description
STDS5D5513	Front Axle Bearings,EXP19LT,1pair, PR
STDS5D5514	Front Axle,EXP19LT,1ea, EA
STDS5D5515	Front Caster,EXP19LT,1ea, EA
STDS5D5516	Front Fork,Blue,EXP19LT,1ea, EA
STDS5D5517	Front Fork,Red,EXP19LT,1ea, EA
STDS5D5518	Fork Stem Bearings,EXP19LT,1/pr, PR
STDS5D1069	Hanger Pins for ATC, SSP, Cruiser,4/pk, BAG
STDS3J0929	Seat Upholstery,19"ATC 3J,5D,9H,1/ea, EA
STDS5D5519	SeatUpholsteryScrews,8Pack,EXP19LT,1pkg, PKG
STDS5D5520	Hand Brake with Cable,EXP19LT,1ea, EA