

Steerable Knee Walker

HCPCS: E0118



- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- New, two-piece padded leg pad provides comfort
- Deluxe braking system
- Removable front basket
- Durable, steel frame with silver vein finish
- Tool free height adjustment
- 7.5" casters are ideal for indoor/outdoor use

PRODUCTS	Description
796	Steerable Knee Walker,Split KneePad,1/EA

SPECIFICATIONS	Description
Warranty	Limited Lifetime, Lifetime Warranty on Brake Cable
Handle Height	31.5"-40.25"
Weight Capacity	350 lb
Carton Shipping Weight	27 lb
Knee Pad Dimensions	7" (W) x 15" (D)
Knee Pad Height	16.75"-20.5"
Weight	21 lb
Wheel Size	7.5"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

PARTS	Description
10230SN	Blow Molded Seat,10230,1ea, EA
9503H79014	Knee Pad,790,1ea, EA
95012R79625	Allen Bolts, 796/RTL799, 4/PKG, PKG
9503H79021	Wheel Bearings,790,PhoenixHD3,1/pr, PR
10230SF	Seat Frame,10230/10233,1/ea, EA
95012R79601	ScrewsForKneePadReplacement, 796, 4/pkg, PKG
95012R79602	Brake Cable Only, 796/RTL799, 1/ea, EA
95012R79603	Wheel Lock Without Pad, Left, 796, 1/ea, EA
95012R79604	Wheel Lock Without Pad, Right, 796, 1/e, EA
95012R79605	Wheel Lock Spring Only, 796, 1/ea, EA

PARTS

Description

95012R79606	BrakePadConnectorBarW/Adj.Bolt, 796,1/e, EA
95012R79607	Tiller Adjustment Knob, 796, 1/ea, EA
95012R79608	Axle Bolt, 796/RTL799, 1/ea, EA
95012R79609	Front/RearWheelW/Bearings&Spacer,796,1e, EA
95012R79610	Handbrake W/ Cable, 796, 1/ea, EA
95012R79612	Knee Pad Only W/Hardware, 796, 1/ea, EA
95012R79613	SeatPostRelLeverW/Hrdwr,796/RTL799, 1/ea, EA
95012R79614	KneePost Adj Pin, 796/RTL799, 1/EA, EA
95012R79615	Brake Pad W/Screws & Nuts, 796, 1/pr, PR
95012R79616	ForkStem Bearing Assy, 796/RTL799, 1/EA, EA