

## Bariatric Sentra EC Heavy-Duty, Extra-Extra-Wide Wheelchair

HCPCS: K0007-1,K0007-1



- Carbon steel frame with triple-coated chrome for an attractive, chip-proof, maintainable finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable reinforced nylon upholstery
- Steel, spoke wheels are strong and maintenance free
- Heavy-duty, flat-free 8" front casters are adjustable in two positions
- Extra-heavy-duty front forks
- Reinforced side frame, caster journals and dual cross bars
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Extra-large, padded armrests provide comfort and support
- Riggings have extra-large, heavy-duty, aluminum foot plates that provide maximum support, are attractive and crack proof (Must be ordered separately, PH-SF for swing-away footrests, PHELR for elevating legrests)
- Comes with push-to-lock wheel locks
- Weight Capacity 700 lbs

### PRODUCTS

### Description

STD26ECDDA	Sentra EC,26"Det.DeskArm,1/cs
STD26ECDFA	Sentra EC,26in Det.FullArm,1/cs
STD28ECDDA	Sentra EC,28"Det.DeskArm,1/cs
STD28ECDFA	Sentra EC,28"Det.FullArm,1/cs
STD30ECDDA	Sentra EC,30"Det.DeskArm,1/cs
STD30ECDFA	Sentra EC,30"Det.FullArm,1/cs

### SPECIFICATIONS

### Description

Warranty	Limited Lifetime
Product Weight	62 lb   64 lb   66 lb
Product Dimensions	43" (L) w/ riggings x 34.5" (W) x 35" (H)
Seat Dimensions	26" (W) x 20" (D)
Weight Capacity	700 lb
Seat To Floor Height	19.5"
Carton Shipping Weight	68 lb   70 lb   72 lb
Folded Dimensions	43" (L) w/ riggings x 14" (W) x 35" (H)
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> ]