

Bariatric Sentra Extra-Heavy-Duty Wheelchair

HCPCS: K0007,K0007



FIGURE A

- Weight Capacity: 500 lbs
- Carbon Steel frame with silver vein finish is attractive, and easy-to-maintain
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable, heavy-gauge reinforced naugahyde upholstery
- Comes standard with carry pocket on backrest
- Reinforced side frame, caster journals and dual cross bars
- Dual-position front fork
- Dual-position front frame
- Composite Mag-style wheels are lightweight and maintenance free
- Hemi-adjustable caster is optional (Part #STDS1074)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Padded armrests provide added user comfort
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Comes with push-to-lock wheel locks

PRODUCTS	Description
STD20DDA-ELR	SentraHD,20"DetDeskArmElvLeg,1/cs
STD20DDA-SF	SentraHD,20"DetDeskArm,Swg.Ft,1/cs
STD22DDA-ELR	SentraHD,22"DetDeskArmElvLeg,1/cs
STD22DDA-SF	SentraHD,22"DetDeskArm,Swg.Ft,1/cs
STD24DDA-ELR	SentraHD,24"DetDeskArmElvLeg,1/cs
STD24DDA-SF	SentraHD,24"DetDeskArm,Swg.Ft,1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	54 lb 56 lb 66.25 lb
Product Dimensions	42" (L) w/ riggings x 29" 31" 32.5" (W) x 37" (H)
Seat Dimensions	20" 22" 24" (W) x 18" (D)
Weight Capacity	500 lb
Seat To Floor Height	19" - 21"
Carton Shipping Weight	60 lb 62 lb 72 lb
Folded Dimensions	42" (L) w/ riggings x 13" (W) x 37" (H)
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov]

ACCESSORIES	Description
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS831	Heel Loops, 1pr, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS803	Oxygn Tank Holdr for W/C,1/cs, CS
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
STDS802	Anti-Tipprs,w/Wheels,Adj,1pr/bx, BX
PARTS	Description
STDS1051	Dust Cap, fork stem, 1/ea, EA
STDS871	Hand Grip,Sentra,1/ea, EA
STDS1026	Front Bumper for Ft.Rst/LegRest,1/ea, EA
STDS1002-1	Front Axle for 8" caster, 1/ea, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
STDS8M4601	ForkStemBearing,SSP2,Cruiser,1/ea, EA
STDS1070	Seat Rail Guide, Removable Arm, 1/ea, EA
STDS4020	Tipping Sleeve 1/ea, EA
STDS005200	SlideTube,Sentra,1/ea, EA
STDS002001	Arm Pad Screw,1/ea, EA
STDS4011	BrakeTip, 1/ea, EA
STDS2329BLK	BackUphol black,SentraHD 20inx 18in,1/bx, BX
STDS2331BLK	Back Uphol 22"x18",Black,1/bx, BX
STDS1052	Rear AxleBolt,SentraHD,1/ea, EA
STDS2333BLK	Back Uphol 24"x18", Black,1/bx, BX
STDS2A5222	Crossframe,22"SentraHD/1ea, EA
STDS1074	6"Hemi Caster,SentraHD,1ea, EA
STDS4010R	Brake Asmb.Right Push to Lock,1/ea, EA
STDS1000HD	Wheel,24""Hvy.Dty.Hand Rl.Sentra,1/ea, EA