



## SHOWER ACCESSORIES

# GRAB BARS

### Features

- Heavy duty steel construction with rust resistant finishes
- Available in chrome plated or white enamel finish
- Knurling provides a good gripping surface
- Mounts vertically, horizontally or diagonally 1 1/2" from wall
- Attractive retail packaging
- Designed to meet the requirements of HCPCS code: E0241
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Three (3) year limited warranty



		Length		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	lb	Kg	
2012A	Chrome Plated Finish, Knurled	12	30.48	250	113.40	1 ea
2016A	Chrome Plated Finish, Knurled	16	40.64	250	113.40	1 ea
2018A	Chrome Plated Finish, Knurled	18	45.72	250	113.40	1 ea
2019A	Chrome Plated Finish, Knurled	24	60.96	250	113.40	1 ea
2021A	Chrome Plated Finish, Knurled	32	81.28	250	113.40	1 ea
40033A	White Versaguard Finish	12	30.48	250	113.40	3/cs
40133A	White Versaguard Finish	18	45.72	250	113.40	3/cs
40233A	White Versaguard Finish	24	60.96	250	113.40	3/c
40333A	White Versaguard Finish	32	81.28	250	113.40	3/cs
40433A	White Versaguard Finish	16	40.64	250	113.40	3/cs



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.