



### Standard Colors:



Red



Blue

## BuzzAbout GP164



Enjoy full-time luxury on the go in the new BuzzAbout! Perfect for travel, the BuzzAbout power wheelchair is in a class all its own in our popular BuzzAround series of mobility products. The BuzzAbout has a stylish captain's seat with contoured support and higher back for longer seating tolerance and comfort. It features motion sensing LED running lights on the footplate and rear handle activated when the joystick is moved. It also offers convenient offboard charging and easy disassembly for transport. The BuzzAbout offers excellent performance and maneuverability indoors and out for your active lifestyle.

### Features:

- Motion Sensing LED Running Lights on Footplate & Rear Handle
- Quick & Easy Disassembly into Four Pieces
- Offboard Charging Capability
- 300 lb. Weight Capacity



### Item

|                                    |               |
|------------------------------------|---------------|
| Description                        | GP164         |
| Medicare Code                      | Not Coded     |
| Weight Capacity                    | 300 lbs.      |
| Drive Wheels                       | Rear          |
| Maximum Speed                      | 4 mph         |
| Operating Range up to <sup>1</sup> | 11 miles      |
| Under Deck Clearance               | 2.5"          |
| Turning Radius                     | 28.5"         |
| Type Batteries                     | 2-22AH        |
| Colors                             | Red, Blue     |
| Freewheel Mode                     | Yes           |
| Electronic Speed Control           | Yes           |
| Electro-Mechanical Brakes          | Yes           |
| Charger                            | Off-board     |
| Controller                         | 40 Amp Remote |

### Dimensions

|  |                      |
|--|----------------------|
| Length                                 | 36"                  |
| Width                                  | 23" at front casters |
| Height (ground to top of back of seat) | 36.75-38.75"         |
| Ground to Top of Seat                  | 21.5"-23.5"          |
| Number of Seat Height Positions        | 3                    |
| Size of Increments                     | 1"                   |

### Weight of Unit: (assembled)

|             |         |
|-------------|---------|
| Front Half  | 29 lbs. |
| Rear Half   | 33 lbs. |
| Seat w/Arms | 34 lbs. |
| Batteries   | 28 lbs. |

### Tires:

|                      |     |
|----------------------|-----|
| Front                | N/A |
| Rear                 | 9"  |
| Casters              | 6"  |
| Rear Anti-Tip Wheels | 2"  |

### Standard Seat

|                           |             |
|---------------------------|-------------|
| Slide Seat                | N/A         |
| Back Height (no headrest) | 18"         |
| Width x Depth             | 17" x 16"   |
| Color                     | Black Vinyl |

### Footrest:

|                   |     |
|-------------------|-----|
| Height adjustable | No  |
| Angle adjustable  | Yes |

### Warranty:

|                          |                   |
|--------------------------|-------------------|
| Frame                    | 5 Years           |
| Drive Train              | 13 Month Warranty |
| Electronics <sup>2</sup> | 13 Month Warranty |
| Batteries                | 1 Year Warranty   |



Front Motion Sensing LED Light



Rear Motion Sensing LED Light



GP164 Shown Disassembled

<sup>1</sup>. Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

<sup>2</sup>. Electronics warranty excludes batteries.