

# So Lite™ S1

## SELL SHEET



### Benefits

- Gives user more freedom to travel with a light weight unit
- Allows user to carry more items with them with three storage options. Two carry bags are removeable and can be taken with you
- Added safety and sense of confidence with electronic stability control feature
- Very comfortable to sit in and intuitive controls

### Product Specifications

Model Reference	08474
Speed (Max)	5 mph (8 km/h)
Battery	1 x 24V 10Ah Lithium
User Weight (Max)	265 lbs (120kg)
Total Weight w/o Batteries	40.8 lbs (18.5kg)
Total Weight with Batteries	44.5 lbs (20.2kg)
Overall Length	1000 mm (39.5 in.)
Overall Width	550 mm (21.5 in.)
Overall Height	940 mm (37 in.)
Folded Length	1000 mm (39 in.)
Folded Width	550 mm (21 in.)
Folded Height	430 mm (17 in.)
Climbing Gradient (Max)	8 Degrees
Turning Radius (Min)	1150 mm (45 in.)
Ground Clearance (Max Mid Frame)	105 mm (4 in.)
Seat Type	Folding with removable seat covers
Seat Width	430 mm (17 in.)
Seat Height	370 mm (18 in.)
Front Tire	150 mm (6 in.) Solid Puncture
Rear Tire	250 mm (10 in.)
Lights	Front LED
Basket	Metal under seat
Tiller Type	Straight Handle
Charger	2 Amp
Charging Type	Battery Box that can be off-board
Brake	Automatic Electromagnetic

### Features

- Very easy tool free assembly
- More storage than other scooter on the market. Comes standard with front and rear carry bag and under seat basket
- Electronic stability control for added safety
- Variety of three colors to choose from
- Holds up to 265 lbs
- Light weight at only 40 lbs (without battery)
- Thick padded seat for extra comfort
- Folds in seconds for easy storage and travel
- Able to be charged by plugging directly into the steering column. Battery removal not required
- Powerful light on the front of the unit