



EVOLUTION
LC-435L



**Storage
Compartment**

*Only available on small,
medium, and large sizes*



**Cup
Holder**

*Only available on small,
medium, and large sizes*

Sitting Pretty™

Evolution Collection



PRIDEMOBILITY.COM

Evolution Collection lift chairs, excluding LC-435 XL and XXL, are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

While Pride Mobility® began exclusively as a lift chair manufacturer in 1986, our company has never stopped innovating, which makes the **Evolution Collection** a fitting name. This collection fits loungers of all kinds, with sizes ranging from small to XXL.

The small, medium, and large sizes feature a cup holder and storage compartment in the pull-out arms, with USB-A and USB-C chargers in the storage area to keep you connected. These sizes also offer an enhanced three-position experience, which provides additional resting positions.

Evolution Collection



Model	LC-435S	LC-435M	LC-435L	LC-435XL	LC-435XXL
Size	Small	Medium	Large	X-Large	XX-Large
User Height	5'3" & under	5'4" - 5'9"	5'9" - 6'	5'8" - 6'1"	5'8" - 6'1"
Weight Capacity	400 lbs.	400 lbs.	400 lbs.	500 lbs.	600 lbs.
Seat-to-Floor	18"	20"	22"	20"	20"
Seat Depth	18"	20"	22"	22"	22"
Seat Width	20"	22"	24"	26"	30"
Top of Back to Seat	24"	26"	28"	29"	29"
Fully Padded Chaise	Yes	Yes	Yes	Yes	Yes
Seating Type	Coil Pocketed Spring	Coil Pocketed Spring	Coil Pocketed Spring	Total Comfort	Total Comfort
Distance from Wall	16"	16.5"	19"	19"	19"
Overall Width	38.5"	40.5"	42"	39.5"	43"
Total Weight	136 lbs.	141 lbs.	146 lbs.	187 lbs.	196 lbs.
Warranty	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated	7-year Pro-rated
Made in America of foreign and domestic parts	No	No	No	Yes	Yes
Position Type	Enhanced 3-position	Enhanced 3-position	Enhanced 3-position	Traditional 3-position	Traditional 3-position

LC-435S / LC-435M / LC-435L FABRIC OPTIONS

Cosi Fabrics




Mocha

Gravel

Blue
Diamond

Syrah

Sta-Kleen¹



Chestnut

LC-435XL / LC-435XXL FABRIC OPTIONS

Cosi Fabrics



Mocha

Gravel

Sta-Kleen¹



Chestnut

Evolution Collection lift chairs, excluding LC-435 XL and XXL, are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

FEATURES

- Easy-to-use toggle hand control
- Three-position lift recliners
 - *Enhanced three-position on small, medium, and large for additional resting positions*
- Integrated battery backup in case of power failure
- Chaise pad design fully supports the legs to eliminate pressure points
- Small, medium, and large models feature a cup holder, armrest storage with USB-A/USB-C charging ports, and magnetic hand control dock on both sides
- XL and XXL models feature a USB charger in the hand control, 2-amp lithium battery backup, and are made in the USA of foreign and domestic parts

OPTIONS & ACCESSORIES


- Lithium battery backup (standard on XL and XXL)
- Head and arm covers
- Non-skid leg levelers
- Cleaning kit

Sitting Pretty™



pridemobility.com

800.800.8586



NOTE: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1) Additional charges will apply for this fabric.