

Transport Chairs

Transport chairs make it easier for caregivers to help individuals get around comfortably and easily. Compared to standard wheelchairs, a transport chair requires a caregiver to push the chair.

A transport chair is much lighter than a wheelchair and comes in a variety of colors and materials such as steel and aluminum. They are more compact, great for light use and getting from one place to another. The wheel sizes vary from 8" to 12" for durability and a 17" or 19" seat width fits everyone's needs.



MDS808200B

Basic Steel Transport Chairs with 8" Wheels

- Durable steel frame with Black upholstery
- Includes detachable swing away footrests for maximum comfort and ease of mobility
- Back folds down for convenient transporting

Item No.	Description	Wt. Capacity	Leg Style	Arm Style	Seat Width	HCPCS	Retail Pkg.	Pkg.
MDS808200B	Black	250 lb (113.4 kg)	Swing-away	Full-Length, Permanent	19"	E0138	N	1 ea
MDS808200BN	Black	300 lb (136.1 kg)	Swing-away	Full-Length, Permanent	17"	–	N	1 ea
MDS808200WAT*	Red	300 lb (136.1 kg)	Swing-away	Desk-Length, Permanent	19"	–	N	1 ea
MDS808200	Chrome	300 lb (136.1 kg)	Swing-away	Full-Length, Permanent	19"	E0138	N	1 ea
New MDS808200BP	Black	250 lb (113.4 kg)	Swing Away	Desk Length, Swing Back	19"	–	N	1 ea
New MDS808200BNP	Black	250 lb (113.4 kg)	Swing Away	Desk Length, Swing Back	17"	–	N	1 ea

*Includes anti-tippers

Aluminum Transport Chairs with 8" Wheels

- Powder-coated aluminum frame with Black upholstery
- Includes detachable swing away footrests and padded full length arm rests

Item No.	Description	Wt. Capacity	Leg Style	Arm Style	Seat Width	HCPCS	Retail Pkg.	Pkg.
MDS808200ABE	Blue	300 lb (136.1 kg)	Swing-away	Full-Length, Permanent	19"	E1038	N	1 ea
MDS808200ARE	Red	300 lb (136.1 kg)	Swing-away	Full-Length, Permanent	19"	E1038	N	1 ea
New MDS808100ARE	Red	300 lb (136.1 kg)	Swing Away	Full Length	17"	–	Y	1 ea
New MDS808100ABE	Blue	300 lb (136.1 kg)	Swing Away	Full Length	17"	–	N	1 ea

Basic Steel Transport Chairs with 12" Rear Wheels

- Durable steel frame with Black upholstery
- Includes detachable footrests and a seat belt for safety
- Back folds down for convenient transporting
- 12" rear wheel for ease on all terrains

Item No.	Description	Wt. Capacity	Leg Style	Arm Style	Seat Width	HCPCS	Retail Pkg.	Pkg.
MDS808210KDR*	Burgundy	300 lb (136.1 kg)	Swing-away	Desk-Length, Permanent	19"	–	N	1 ea
MDS808200KD	Red	300 lb (136.1 kg)	Swing-away	Desk-Length, Permanent	19"	–	Y	1 ea
MDS808200KDT*	Teal	300 lb (136.1 kg)	Swing-away	Full-Length, Permanent	19"	–	Y	1 ea
New MDS808200KWM	Tan	300 lb (136.1 kg)	Swing Away	Desk Length	18.5"	–	N	1 ea
New MDS808200KNM	Navy	300 lb (136.1 kg)	Swing Away	Desk Length, Swing Back	17"	–	N	1 ea

*Antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

Aluminum Transport Chairs with 12" Wheels

- Powder coated aluminum frame with Black upholstery
- Includes detachable swing away footrests and padded full length arm rests
- Larger back tire for maneuvering through obstacles
- Hand brakes on handles

Item No.	Description	Wt. Capacity	Leg Style	Arm Style	Seat Width	HCPCS	Retail Pkg.	Pkg.
MDS808210ABE	Blue	300 lb (136.1 kg)	Swing-Away	Full-Length, Permanent	19"	E1038	N	1 ea
MDS808210ARE	Red	300 lb (136.1 kg)	Swing-Away	Full-Length, Permanent	19"	E1038	N	1 ea