

# VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST  
COMFORTABLE  
RECLINER  
IN THE WORLD!

## THE 4-MOTOR RECLINER

NOW  
AVAILABLE

TRUE-INFINITE  
POSITIONING AND  
MEMORY TOGGLE  
HAND CONTROL!

FDA Class II Medical Device¹

Shown in Astro Brown Fabric



COUNTLESS  
POSITION  
COMBINATIONS

## TRANQUIL 2¹

The **Tranquil 2** combines quiet comfort with gorgeous fabrics, for a rich and upscale look. Enhanced recline positioning to raise your feet above your heart.

The Tranquil 2 has been upgraded with an ergonomically designed memory toggle remote that features two programmable memory position buttons along with up/down/recline functions and custom comfort levels.



Scan here for  
more information  
on **VivaLift!®**









**PRIDE**®  
MOBILITY  
LIVE YOUR BEST®

# TRANQUIL 2<sup>1</sup>

**VivaLIFT!**<sup>®</sup>  
POWER RECLINER COLLECTION<sup>1</sup>

COLLECTION	TRANQUIL 2 <sup>1</sup>	TRANQUIL 2 <sup>1</sup>	TRANQUIL 2 <sup>1</sup>
MODEL	PLR-935S	PLR-935M	PLR-935LT
POSITION	True-Infinite Position	True-Infinite Position	True-Infinite Position
HEIGHT RANGE	5'3" and below	5'4"-5'9"	5'10" and above
SIZE	Small	Medium	Large/Tall
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	18"	19.5"	21.5"
SEAT DEPTH	18.5"	20.5"	22.5"
SEAT WIDTH	20"	22"	23.5"
TOP OF BACK TO SEAT	25"	27"	30"
FULLY PADDED CHAISE	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	40"	42"	47"
DISTANCE FROM WALL	13"	16"	16"
OVERALL WIDTH	33.75"	37.5"	38"
CHAIR WEIGHT	137 lbs.	145 lbs.	155 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

## INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. **NEW!** Enhanced true-infinite position recline enables the backrest and footrest to move independently so that a variety of positions can be reached including raising your feet above your heart.

## STANDARD FEATURES

- 7. NEW!** Ergonomically designed toggle remote has two programmable memory position buttons. Performs the up/down, recline functions. Additionally, it has a convenient USB charger port that allows you to charge your smart device.
- 8.** Lithium battery backup ensures your recliner is operational when the power goes out.

1. VivaLift!<sup>®</sup> Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



**ES** = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

## STANDARD FABRIC: ASTRO



GREY

BROWN

MUSHROOM

(US) 800-800-8586 / (CAN) 888-570-1113 / [pridemobility.com](http://pridemobility.com)

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>