





Standard Fabrics:



Acorn



Graphite



Maple

This faux leather is a rich and modern alternative to natural leather that is extremely comfortable while offering the rugged durability of leather. Made from 100% polyester, the fabric is easy to care for ensuring functionality as well as elegance. This long-lasting fabric will complement a variety of decors while standing the test of time under heavy use.

Crossover fabrics are unavailable.

Dione PR446



Comfort Zone™
by GOLDEN

The Dione is offered in both medium and large sizes and features a stylish, overstuffed seamed back and plush chaise pad. It comes standard with adjustable lumbar and headrest, programmable AutoDrive™ hand control, and USB charging port for your phone or tablet. Also standard is our patented four Comfort Zone positioning system, which provides a stress-free and supportive reclining experience.



Standard Features:

- Patented MaxiComfort® Positioning
- Adjustable Headrest
- Adjustable Lumbar
- Programmable AutoDrive™ Hand Control
- USB Charging Port
- Extra Pocket
- SmartTek® Diagnostic System



DeLUNA Series PR446

Model Number	PR446-MED	PR446-LAR
Size	Medium	Large
Class II Medical Device	YES	YES
Weight Capacity	375 LBS	375 LBS
Positioning	MaxiComfort	MaxiComfort
Comfort Zone(s)	4	4
Twilight Technology	NO	NO
ZG+	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	YES	YES
Back Type	Seam	Seam
Removable Back Type	Quick Release	Quick Release
Head Rest Cover	NO	NO
Arm Rest Covers	NO	NO
Overall Width	34"	35"
Overall Height	43.5"	46"
Seat Width	22"	22.5"
Floor to Top of Seat	20.5"	22"
Seat Depth	21"	21.5"
Seat to Top of Back	28"	30"
Seat to Top of Arm	7"	7.5"
Distance Required from Wall Reclined	24"	26"
Extension Reclined from Back of Seat to End of Foot Rest	42.5"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Arm	Inside Arm
Extra Pocket Option	Standard	Standard
Heat & Vibrating Massage Option	NO	NO
HeatWave Technology Option	NO	NO
Articulating Headrest Option	Standard	Standard
Power Adjustable Lumbar Option	Standard	Standard
Foot Rest Extension Option	NO	NO
Straight Lift Option	Standard	Standard
USB Charging Port	Right Arm	Right Arm
Battery Backup	Lithium-ion	Lithium-ion
SmartTek Diagnostic System	Standard	Standard
Place of Origin	Imported	Imported

DeLuna™
SERIES



Shown
in Graphite



Convenient built-in USB charging port for your phone or tablet located on the outside arm.

Designed and engineered in Old Forge, PA. Assembled in China using Golden components.



Comfort Zone™
by GOLDEN

Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.

GOLDEN

Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68

Find us on
Facebook

www.facebook.com/goldentechnologies

SHT-PR446-010924