

|                                                            |                                                                                                                                                          |
|------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Model</b>                                               | Maxima <sup>1</sup>                                                                                                                                      |
| <b>Part Number</b>                                         | 3W: SC901 4W: SC941                                                                                                                                      |
| <b>Weight Capacity</b>                                     | 500 lbs.                                                                                                                                                 |
| <b>Maximum Speed<sup>2</sup></b>                           | Up to 5.8 mph                                                                                                                                            |
| <b>Ground Clearance<sup>3</sup></b>                        | 3W: 2.5" at motor; 4W: 2.25" at motor                                                                                                                    |
| <b>Turning Radius<sup>3</sup></b>                          | 3W: 49"; 4W: 53"                                                                                                                                         |
| <b>Overall Length<sup>3</sup></b>                          | 47.5"                                                                                                                                                    |
| <b>Overall Width<sup>3</sup></b>                           | 27" (not including armrest)                                                                                                                              |
| <b>Seat-to-Ground Height<sup>3</sup></b>                   | 3W: 22.76" – 24.76";<br>4W: 22.76" – 24.76"<br>w/ PES: 21" – 27.25"                                                                                      |
| <b>Seat-to-Deck Height<sup>3</sup></b>                     | 3W: 16.25"–18.25"; 4W: 16"–18"<br>w/ PES: 16"– 22.25"                                                                                                    |
| <b>Front Tires<sup>3</sup></b>                             | 3W: 10.5" x 3.6" solid<br>4W: 10" x 3.6" solid                                                                                                           |
| <b>Rear Tires<sup>3</sup></b>                              | 10.75" x 3.6" solid                                                                                                                                      |
| <b>Suspension</b>                                          | N/A                                                                                                                                                      |
| <b>Range Per Charge<sup>2,5</sup></b>                      | 3W: Up to 19.8 miles<br>4W: Up to 19.5 miles                                                                                                             |
| <b>Total Weight With and Without Batteries<sup>7</sup></b> | 3W (w/out batteries): 183.5 lbs.<br>3W (w/batteries): 256.5 lbs.<br>4W (w/out batteries): 196.8 lbs.<br>4W (w/batteries): 269.8 lbs.                     |
| <b>Heaviest Piece When Disassembled<sup>3</sup></b>        | N/A (unit does not disassemble)                                                                                                                          |
| <b>Standard Seating<sup>3</sup></b>                        | Type: Depth-adjustable, high-back<br>Weight: 53.2 lbs.<br>Material: Black vinyl<br>Dimensions: 22"x 18"-20"<br>(depth adjustable)<br>Back height: 19.75" |
| <b>Drive System</b>                                        | Sealed transaxle, 24-volt DC motor                                                                                                                       |
| <b>Dual Braking System</b>                                 | Regenerative and electromechanical                                                                                                                       |
| <b>Battery Requirements<sup>4,6</sup></b>                  | Size: (2) Gel, 12V, 58Ah, NF-22<br>Weight: 36.5 lbs. each                                                                                                |
| <b>Battery Charger</b>                                     | Off-board, 5-amp                                                                                                                                         |
| <b>Warranty</b>                                            | Frame: lifetime limited<br>Drivetrain: 2-year limited<br>Electronics: 2-year limited<br>Batteries: 6-month limited                                       |
| <b>Color Option</b>                                        | Garnet Red (matte)<br>Ocean Blue (matte)                                                                                                                 |

3W: HCPCS Code: K0802 (Group 1, Very Heavy Duty)  
4W: HCPCS Code: K0808 (Group 2, Very Heavy Duty)

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.
- 7 Includes standard seat weight.

**PRIDE  
MOBILITY®**  
800-800-8586 • [pridemobility.com](http://pridemobility.com)



HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

092525

# MAXIMA<sup>1</sup>

## Ultra Heavy-Duty for Maximum Performance



FDA Class II  
Medical Device<sup>1</sup>

800-800-8586 • [pridemobility.com](http://pridemobility.com)



# HEAVY-DUTY PERFORMANCE



## FEATURES

Infinite position tiller lever (mech lock)

22" x 18"-20" High back depth adjustable seat

- 106° limited recline
- 4-point swivel (no sliders)
- 7" headrest with 3" of adjustment

USB charger built into the tiller for convenient charging of smart phones and portable electronic devices

Full LED light package includes front headlights, turn signals and rear taillights

Ambient voltmeter screen brightness automatically adjusts based on light intensity

Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand

Longer per-charge range

- Standard with front and side bumper (3W) and side bumper (4W)
- XLR battery charger on tiller

## CONSOLE

- Ambient voltmeter
- Horn
- Light switch
- On/off key switch
- Speed control
- Turn signal button

## OPTION

- Power elevating seat available on 3 & 4 wheel
- 6.25" of additional height to access cabinets, countertops and other out-of-reach items
- High back, black vinyl 20" x 20" seat w/o sliders

## AVAILABLE COLORS


Garnet Red (matte)

Ocean Blue (matte)

## OPTIONAL ACCESSORIES\*

- Additional rearview mirror
- Cup holder (B)
- Front bumper (4w)
- Oxygen holder
- Rear basket
- Rear bumper (3w, 4w)
- Saddle bag
- Safety flag
- Single/double cane/crutch holder
- Weather covers
- Military Patches

\* Only one rear accessory can be used at a time.



Full LED light package



Infinite position tiller lever (mech lock)



3W — FDA Class II Medical Device<sup>1</sup>