

Model	Victory® 10.2 ¹
Product Number	SC6102
Weight Capacity	400 lbs.
Maximum Speed ²	Up to 5.2 mph
Ground Clearance ⁴	2.4" at transaxle, 2.1" at anti-tip 3" at side deck
Turning Radius ⁴	45"
Overall Length ⁴	45"
Overall Width ⁴	22.25"
Seat-to-Ground Height	22" - 24"
Seat-to-Deck Height	15.25"-17.25"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	N/A
Range Per Charge ^{2,3}	40 AH: Up to 13 miles U-1: Up to 11 miles
Total Weight ⁶	Without batteries: 130.25 lbs. With batteries: 178.25 lbs. (U1); 194.25 lbs. (40AH)
Heaviest Piece When Disassembled ⁴	56.5 lbs. (front section)
Standard Seating	Type: Foldable Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width, 17" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements ^{5,7}	Size: (2) U1 standard, Weight: 24.5 lbs. each Size: (2) 40 AH (optional), Weight: 32.5 lbs. each
Battery Charger	3.5-amp, off-board (U1), 5-amp, off-board (40 AH)
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

PRIDE MOBILITY®

800-800-8586 • pridemobility.com



HCPCS Code: K0801 (Group 1, Heavy Duty)

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ AGM or gel cell type required.
- ⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁶ Battery weight may vary based on manufacturer.
- ⁷ Includes standard seat weight.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

082525

VICTORY® 10.2

CONTEMPORARY FEATURES
ENHANCE YOUR ACTIVE LIFESTYLE



FDA Class II Medical Device¹

800-800-8586 • pridemobility.com



CONVENIENCE AND VERSATILITY



Available in
Candy Apple Red

CONVENIENCE AND VERSATILITY

The Victory® 10.2 offers great convenience and versatility with exclusive feather-touch disassembly and a delta tiller with wrap-around handles. Contemporary features such as the integrated storage/cup holder on the tiller ensure ample storage options, while a new under-tiller puddle light offers improved visibility to both the storage area and USB charging port.

FEATURES

- Exclusive one-hand feather touch disassembly
- Integrated tiller charger port lets you conveniently charge your scooter and mobile devices
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Integrated storage/cup holder in tiller

Optional Seating

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests:

- 18" x 18" (w/sliders)
- 20" x 20" • 22" x 18"

Deluxe contour high back black vinyl with height adjustable headrest, flip-up height and width adjustable armrests (no sliders):

- 18" x 18" • 18" x 20"
- 20" x 20" • 22" x 20"

AVAILABLE COLORS

Candy Apple Red

FDA Class II Medical Device¹



Frosted LED headlights



Under-tiller puddle light



Ambient voltmeter

OPTIONAL ACCESSORIES*

- 50", 60", 70" Lap belt
- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Military Patches
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover

* Only one rear accessory can be used at a time.