

# ROYAL ALUMINUM ROYAL STEEL

## 300 lb Weight Capacity



ITEM	COLORS
807S - Royal Steel	807SBB (Blueberry) 807SGT (Garnet)
807 - Royal Aluminum	807BB (Blueberry) 807RD (Red) 807SG (Slate Grey)

SPECIFICATIONS	807	807S
Weight	14 lbs	17 lbs
Handle Height	30" - 37"	32" - 38"
Height Adjustment	Handles only	Handles only
Seat to Floor Height	20.25"	22"
Seat Dimensions	17.5" W x 13" D	17.5" W x 13" D
Unit Dimensions	23.25" W x 26" L	23.25" W x 26" L
Frame	Aluminum Round	Steel Round
Wheels	6"	6"
Weight Capacity	300 lbs	300 lbs

## WARRANTY

Rhythm Healthcare warrants this product to be free of defects in materials and workmanship for the lifetime of the original consumer. We offer this **Limited Lifetime Warranty** to affirm our confidence that the materials and workmanship of our products will result in years of dependable service.

This warranty is conditioned upon the proper use of this product in accordance with the instructions and manuals, including proper use and maintenance. This warranty does not cover non-durable components, such as rubber accessories, which are subject to normal wear and will require periodic replacement. In the event warranty service is required, we will determine whether the product is to be repaired or replaced.

For questions or to obtain warranty service, please contact an authorized Rhythm Healthcare dealer.