



VGuard®

Vinyl Exam Gloves




Features and Benefits

- Economical, all purpose glove.
- Clear Color.
- Powder Free.
- Medical Grade.
- Single Use Only.
- Smooth Surface.
- Prop 65 Compliant.
- Complies with FFDCA food handling requirements 21 CFR, Parts 170-199.

Specifications



Usage



Ordering Information

ITEM NUMBER	SIZE	FINGER THICKNESS	PALM THICKNESS	GRAM WEIGHT*	INNER BOX	OUTER CASE
A21A11	SM	4.0 mil +/- 0.5 mil	3.0 mil +/- 0.5 mil	5.0g +/- 0.2g	8-50001-98101-9	8-50001-98150-7
A21A12	MD				8-50001-98103-3	8-50001-98104-0
A21A13	LG				8-50001-98106-4	8-50001-98107-1
A21A14	XL				8-50001-98109-5	8-50001-98110-1

*Based on size M

ITEM NUMBER	QUANTITY	INNER BOX DIMENSIONS	CASE DIMENSIONS	TI/HI	CASE WEIGHT (LBS.)	COUNTRY OF ORIGIN
A21A11	100/bx 10 bx/cs	9.1 x 4.7 x 2.6"	13.4 x 9.8 x 9.6"	14/6	11.9	China
A21A12					13.0	
A21A13					14.1	
A21A14					15.2	

Finding the right hand protection doesn't have to be complicated. Consider VGuard® gloves as an extension of our promise to you: High quality products, great value and excellence in service. If ever you have a question regarding which glove is best for your application, we're here to help. VGuard® is protection you can count on!

Vanguard Safety

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)
Info@VanguardSafety.com • www.VanguardSafety.com

© 2024 VGuard® is a registered trademark of Vanguard Safety



Vinyl Exam Gloves



Composition

Material:	Polyvinyl Chloride (OVS)
Allergen content:	Contains no proteins found in latex
Powdered:	No, as per ASTM D6124

Quality

Quality Control:	ASTM D5250, ASTM F1671
Quality Assurance:	ISO:9001, ISO 13485, GMP, HACCP
Regulatory Compliance:	FDA approved

Technical Data

THICKNESS	MIL
Palm:	3.0 ± 0.5
Finger:	4.0 ± 0.5

Weight*:	5.0 gm ± 0.2 gm
-----------------	-----------------

*Based on size M

PHYSICAL PROPERTIES	BEFORE AGING	AFTER AGING
Tensile Strength (MPa):	Min. 11	Min. 11
Ultimate Elongation (%):	Min. 300	Min. 300
Freedom from holes:	AQL 2.5	N/A

DIMENSIONS			
ITEM NO.	SIZE	PALM WIDTH (mm)	LENGTH (mm)
A21A11	Small	80 ± 10	Min. 220
A21A12	Medium	95 ± 10	Min. 230
A21A13	Large	110 ± 10	Min. 230
A21A14	XL	120 ± 10	Min. 230

COUNTRY OF ORIGIN	HARMONIZATION CODE
China	3926.20.1010

Vanguard Safety

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)
Info@VanguardSafety.com • www.VanguardSafety.com