



VGuard®

Nitrile Chemo Exam Gloves




Specifications



Features and Benefits

- Effective protection against chemicals, oils and bloodborne pathogens.
- Excellent puncture and abrasion resistance.
- Medical Grade. Chemotherapy Rated.
- Powder Free.
- Ambidextrous.
- Eliminates exposure to latex.
- Textured surface for improved grip.
- Twice the gloves per case for reduced packaging waste. Fits standard glove dispensers.

Usage



Ordering Information

ITEM NUMBER	SIZE	FINGER THICKNESS	PALM THICKNESS	GRAM WEIGHT*	INNER BOX	OUTER CASE
A14C11	SM	3.5 mil +/- 0.5 mil	2.8 mil +/- 0.5 mil	3.5g +/- 0.2g	008-50001-98172-9	008-50001-98173-6
A14C12	MD				008-50001-98175-0	008-50001-98176-7
A14C13	LG				008-50001-98178-1	008-50001-98179-8
A14C14	XL				008-50001-98181-1	008-50001-98182-8

*Based on size M

ITEM NUMBER	QUANTITY	INNER BOX DIMENSIONS	CASE DIMENSIONS	TI/HI	CASE WEIGHT (LBS.)	COUNTRY OF ORIGIN
A14C11	200/bx 10 bx/cs	9.7 x 4.7 x 3.2"	16.3 x 10.4 x 10"	10/6	16.8	China Malaysia Thailand
A14C12					18.1	
A14C13					19.4	
A14C14					20.7	

Finding the right hand protection doesn't have to be complicated. Consider VGuard® gloves as an extension of our promise to you: High quality products, great value and excellence in service. If ever you have a question regarding which glove is best for your application, we're here to help. VGuard® is protection you can count on!

Vanguard Safety

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)

Info@VanguardSafety.com • www.VanguardSafety.com

© 2025 VGuard® is a registered trademark of Vanguard Safety



Nitrile Chemo Exam Gloves



Composition

Material:	Nitrile synthetic rubber
Allergen content:	Contains no proteins found in latex
Powdered:	No, as per ASTM D6124

Quality

Quality Control:	ASTM D6319, ASTM F1671 ASTM D6978-05 (Chemo)
Quality Assurance:	ISO: 9001, ISO:14001, ISO:13485, GMP, HACCP
Regulatory Compliance:	FDA approved

Technical Data

THICKNESS	MIL
Palm:	2.8 ± 0.5
Finger:	3.5 ± 0.5

Weight*:	3.5 gm ± 0.2 gm
-----------------	-----------------

*Based on size M

PHYSICAL PROPERTIES	BEFORE AGING	AFTER AGING
Tensile Strength (MPa):	Min. 14	Min. 14
Ultimate Elongation (%):	Min. 500	Min. 400
Freedom from holes:	AQL 2.5	N/A

DIMENSIONS			
ITEM NO.	SIZE	PALM WIDTH (mm)	LENGTH (mm)
A14C11	Small	80 ± 10	Min. 220
A14C12	Medium	95 ± 10	Min. 230
A14C13	Large	110 ± 10	Min. 230
A14C14	XL	120 ± 10	Min. 230

COUNTRY OF ORIGIN	HARMONIZATION CODE
China Malaysia Thailand	4015.12.1010

Vanguard Safety

Cumming, GA 30041 • 1-866-VGUARD-6 (1-866-848-2736)
Info@VanguardSafety.com • www.VanguardSafety.com