

Model	Victory® 10.2 ¹
Product Number	SC6102
Weight Capacity	400 lbs.
Maximum Speed²	Up to 5.2 mph
Ground Clearance⁴	2.4" at transaxle, 2.1" at anti-tip 3" at side deck
Turning Radius⁴	45"
Overall Length⁴	45"
Overall Width⁴	22.25"
Seat-to-Ground Height	22" - 24"
Seat-to-Deck Height	15.25" - 17.25"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	N/A
Range Per Charge^{2,3}	40 AH: Up to 13 miles U-1: Up to 11 miles
Total Weight⁶	Without batteries: 130.25 lbs. With batteries: 178.25 lbs. (U1); 194.25 lbs. (40AH)
Heaviest Piece When Disassembled⁴	56.5 lbs. (front section)
Standard Seating	Type: Foldable Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width, 17" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements^{5,7}	Size: (2) U1 standard, Weight: 24.5 lbs. each Size: (2) 40 AH (optional), Weight: 32.5 lbs. each
Battery Charger	3.5-amp, off-board (U1), 5-amp, off-board (40 AH)
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

PRIDE
MOBILITY

(US) 800.800.8586 • pridemobility.com



HCPCS Code: K0801 (Group 1, Heavy Duty)

- ¹ Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- ² Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ⁴ AGM or gel cell type required.
- ⁵ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁶ Battery weight may vary based on manufacturer.
- ⁷ Includes standard seat weight.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

051324

VICTORY® 10.2

CONTEMPORARY FEATURES
ENHANCE YOUR ACTIVE LIFESTYLE



FDA Class II Medical Device¹

LIVE YOUR BEST®

CONVENIENCE AND VERSATILITY



Available in
Candy Apple Red

18" x 17"
Pinchless hinge,
black vinyl seat
with 8-point
swivel and
sliders for
increased
comfort

Integrated
storage/cup
holder and USB
charging port

Standard
front basket

Easy-to-grip
tiller adjustment
knob

Long-lasting,
maximum-
intensity LED
headlight
offers optimal
pathway
illumination

Exclusive
stylish,
lightweight,
non-scuffing,
black, low-
profile tires

FDA Class II Medical Device¹

CONVENIENCE AND VERSATILITY

The Victory® 10.2 offers great convenience and versatility with exclusive feather-touch disassembly and a delta tiller with wrap-around handles. Contemporary features such as the integrated storage/cup holder on the tiller ensure ample storage options, while a new under-tiller puddle light offers improved visibility to both the storage area and USB charging port.

Features

- Exclusive one-hand feather touch disassembly
- Integrated tiller charger port lets you conveniently charge your scooter and mobile devices
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Integrated storage/cup holder in tiller

Accessories

- 50", 60", 70" Lap belt
- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Military Patches
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder

Optional Seating

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests:

- 18" x 18" (w/sliders)
- 20" x 20" • 22" x 18"

Deluxe contour high back black vinyl with height adjustable headrest, flip-up height and width adjustable armrests (no sliders):

- 18" x 18" • 18" x 20"
- 20" x 20" • 22" x 20"



Frosted LED headlights

Under-tiller puddle light

Ambient voltmeter