

|   |   |
|---|---|
| Model   | Victory® LX Sport <sup>1</sup>  |
| Product Number                                | 4W: S710LXW   |
| Weight Capacity                               | 400 lbs.  |
| Maximum Speed <sup>2</sup>                    | Up to 8.6 mph (200 lbs.)<br>Up to 8 mph (400 lbs.)  |
| Ground Clearance <sup>4</sup>                 | 2.2" at motor, 3.15" at mid-frame   |
| Turning Radius <sup>4</sup>                   | 61.25"  |
| Overall Length <sup>4</sup>                   | 47"   |
| Overall Width <sup>4</sup>                    | 25.5"   |
| Seat-to-Ground Height                         | 23" - 25"   |
| Seat-to-Deck Height                           | 16.5" -18.5"  |
| Front Tires                                   | 3.5" x 10" solid  |
| Rear Tires                                    | 3.5" x 10" solid  |
| Suspension                                    | Front and rear  |
| Range Per Charge <sup>2,3</sup>               | Up to 19.7 miles at 200 lbs.<br>Up to 13 miles at 400 lbs.  |
| Total Weight <sup>6</sup>                     | Without batteries: 160 lbs.<br>With batteries: 220.4 lbs.   |
| Heaviest Piece When Disassembled <sup>4</sup> | 77 lbs. (front section)   |
| Standard Seating                              | Type: Limited reclining high-back<br>Weight: 43 lbs.<br>Material: Black Vinyl<br>Dimensions: 18" width, 18" depth |
| Drive System                                  | Rear-wheel drive, sealed transaxle with 24-volt DC motor  |
| Dual Braking System                           | Regenerative and electromechanical  |
| Battery Requirements <sup>5,7</sup>           | Size: (2) 50 AH<br>Weight: 30.2 lbs. each   |
| Battery Charger                               | 5 amp, off-board  |
| Warranty                                      | Frame: lifetime limited<br>Drivetrain: 2-year limited<br>Electronics: 2-year limited<br>Battery: 6-months         |

## HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

- <sup>1</sup> Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- <sup>2</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- <sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- <sup>4</sup> AGM or gel cell type required.
- <sup>5</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>6</sup> Battery weight may vary based on manufacturer.
- <sup>7</sup> Includes standard seat weight.



FDA Class II Medical Device<sup>1</sup>

**PRIDE**  
MOBILITY  
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113

FOLLOW US        
pridemobility.com

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

101524

**VICTORY<sup>®</sup> LX<sup>1</sup>**  
**SPORT**

ENJOY A  
**SMOOTH**  
AND  
**SPORTY**  
**RIDE**



FDA Class II Medical Device<sup>1</sup>

**LIVE YOUR BEST®**

FOLLOW US        
pridemobility.com

# SMOOTH, COMFORTABLE RIDE

18" x 18"  
Vinyl reclining  
high-back seat  
with sliders,  
height adjustable  
headrest and  
flip-up, height  
and width-  
adjustable  
armrests

Puddle light  
under tiller

Standard  
front basket

Easy-to-grip  
tiller adjustment  
knob

Full LED  
lighting  
package  
includes  
curb light,  
rear lights  
and  
directional  
signals

Storage  
compartment  
under the tiller

CTS Suspension includes adjustable  
shocks for greater comfort



Tiller console



CTS Suspension

Maintain your active lifestyle with the **Victory® LX Sport**. Equipped with **Comfort-Trac Suspension (CTS)**, the Victory LX Sport easily handles more varied terrain while providing a smooth and comfortable ride. Enjoy convenient storage and a courtesy light under the tiller. A long-lasting LED headlight ensures optimal pathway illumination. Experience great features and a sporty style with the Victory LX Sport.

- CTS Suspension includes adjustable shocks for greater comfort
- Exclusive one-hand feather touch disassembly
- Charger port located on the tiller lets you conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Easy-to-read LED voltmeter
- Rearview mirror

## Colors



True Blue

Candy Apple Red

## Optional Seating

Limited recline,  
black vinyl, high  
back with height  
adjustable  
headrest, flip-up  
height and width  
adjustable  
armrests:

- 20"x 20"

Black vinyl,  
high back without  
sliders:

- 22"x18"

## Accessories

- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Lap belts: 50", 60" and 70"

## Military Patches

