



VIVE BACK BRACE

SUP1040

HCPCS: L0625

Ergonomically designed to effectively target lower back pain, the Vive universal back brace supports the lumbar region with adjustable compression, reinforced vertical stays and a removable cushioned lumbar pad. The lightweight, breathable neoprene wrap is nonrestrictive, allowing you to easily move throughout your day. Adjusting to fit most waist sizes, the compression brace attaches with a strong dual fastener system that will securely hold the brace in place without rolling, slipping or bunching.

SPECIFICATIONS

Materials:

- Breathable neoprene.
- Hook and loop fastening material.

Care and cleaning:

- Hand wash cold. No bleach.
- Air dry.

Sizes:

- Standard Size: 26" - 44"
- X-Large Size: 44" - 62"

What's included:

- Lower back brace
- Lumbar Insert
- 60 day guarantee

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



SUPPORT FOR YOUR LOWER BACK

Provides all-day comfort & pain relief.



SECURE, NON-SLIP FIT

For your perfect fit during any activity.



FIT BOTH MEN & WOMEN

Adjustable compression for waist sizes up to 44".

SUPPORT FOR LOWER BACK PAIN RELIEF

Providing targeted support for the lower back, the Vive universal back brace relieves chronic pain and pain caused by sciatica, degenerative disk disease, slipped or herniated disks, scoliosis, spasms and injury. The adjustable compression wrap is reinforced with integrated vertical stays for optimal support without restricting movement.

REMOVABLE LUMBAR PAD

An ergonomically designed lumbar pad supplies additional pressure and support to the lower lumbar region, providing targeted relief from back pain. The cushioned lumbar pad is also removable for easy customization. Reduces the risk of reinjury with gentle compression therapy that retains therapeutic heat to promote healing by increasing circulation.

DUAL ADJUSTABLE FASTENERS

Easily adjust the level of compression produced with the dual fasteners. Constructed with a strong latching material for a secure, nonslip fit, the brace is adjustable, accommodating waist sizes between 26" to 44".

BREATHABLE AND LIGHTWEIGHT

Breathable neoprene brace is lightweight and durable. Naturally moisture wicking, the compression wrap stays comfortably dry all day and is washable for added convenience. The nonslip material will not bunch, pinch or irritate the skin.



(239) 220-5772



sales@vivehealth.com

vive
health

SR23 - 7R1