

VIVE PRECISION DIGITAL THERMOMETER

DMD1055WHT

Experience quick and reliable results with our versatile Vive Precision Digital Thermometer. It offers a one-touch operation for fast temperature readings, along with audible fever alerts for rapid fever detection. The thermometer is rigorously tested for guaranteed accuracy and conveniently recalls the last reading for easy temperature tracking. You can choose between Celsius and Fahrenheit readings and use it orally, rectally, or axillary. Plus, it's easy to clean and sanitize, reducing the risk of cross-contamination. With an out-of-the-box ready design, it comes with a battery and protective storage case for cleanliness and convenience.

SPECIFICATIONS

What's Included:

- ♥ 1 Vive Digital Thermometer
- ♥ One LR41 battery
- ♥ Protective case

Measurement Range:

- ♥ Celsius: 32° - 42.9°
- ♥ Fahrenheit: 90° - 109.9°

Display:

- ♥ Digital LCD display
- ♥ Displays up to four digits
- ♥ Low battery indicator

Power:

- ♥ Requires one 1.5V alkaline battery (LR41)
- ♥ Battery life: One year or 1200 temperature reading

Care: Store thermometer in a dry location free from dust and direct sunlight. Dampen a cloth with a disinfectant (e.g. ethyl alcohol) and wipe clean.





QUICK RELIABLE RESULTS

Temperature readings are quickly displayed with audible fever alerts for fast fever detection.

GUARANTEED ACCURACY

Each thermometer is rigorously tested to ensure clinically accurate readings and recalls the last reading for easy temperature tracking.

VERSATILE THERMOMETER

Displaying temperature readings in Celsius or Fahrenheit, the thermometer can be used orally, rectally, or axillary.

ONE-TOUCH OPERATION

Simple to use with a one-touch temperature reading and battery-saving auto-off function.

EASY TO CLEAN

The digital thermometer is easy to clean and sanitize to reduce cross-contamination.

OUT OF THE BOX READY

Includes the digital thermometer, battery, and protective storage case to keep the thermometer clean and ready to use.



vive now

- ♥ Unlimited data tracking
- ♥ Automatic measurement transfer
- ♥ Exportable PDF summaries and complete CSV file
- ♥ Graphical display of measurements
- ♥ Average values for specific date ranges



CALL NOW!



www.vivenow.health



(239) 228-8939

CALL NOW!

SR24 - BR1