



VIVE **FOREARM CRUTCHES**

MOB1014BLK

HCPCS: E0110

Sleek and stylish, the euro design forearm crutches by Vive offer stabilizing support and greater mobility when recovering from an injury or surgery or managing a long term disability. Ergonomically designed arm cuffs and hand grips provide a comfortable experience throughout the day. Lightweight and durable, the corrosion-resistant aluminum frames are also height adjustable for a customized fit.

SPECIFICATIONS

Materials:

➤ Lightweight aluminum frame. Composite arm cuff and handgrip. Rubber nonslip crutch caps.

Height adjustable:

➤ From 31" to 40" at the handgrip.

Care Instructions:

➤ Wipe down crutches with a damp cloth and mild detergent as needed. Dry immediately with a lint-free cloth.

What's included:

➤ Forearm crutches - one pair.

Color:

➤ Black

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



ERGONOMIC DESIGN

Eliminates underarm pain caused by traditional crutches.



CONTOURED GRIP

Effectively reduces hand and wrist fatigue.



ADJUSTABLE HEIGHT

For comfortable and customized fit.

DON'T LET YOUR INJURY HOLD YOU BACK

Providing a comfortable and secure experience, the Vive Forearm Crutches eliminate the fatigue and pain caused by standard crutches while allowing you independent mobility when recovering from an injury or surgery. The supportive, ergonomic grip is contoured to fit the hand, reducing hand and wrist fatigue and ensuring proper positioning for safe mobility throughout any activity.

STYLISH NEW FINISH

The Vive Forearm Crutches have a stylish steel gray matte finish with black arm cuffs and hand grips. The molded arm cuff and hand grips will not discolor your clothing.

ADJUSTABLE LENGTH FOR COMFORTABLE AND CUSTOMIZED FIT

Easy push button adjustments allow the height to be adjusted in 1" increments for a customized fit. Each crutch is capped with a non-skid, non-marking rubber tip for added safety and security.

ONE DURABLE AND LIGHTWEIGHT PAIR

Arm cuff and hand grip are molded as one piece for greater stability and durability. The lightweight aluminum frame allows you to move easily. Corrosion-resistant, the durable frame supports up to 300 pounds.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - 2RI