



VIVE

METAL GRAB BAR

LVA1079

The Vive grab bar offers a secure, comfortable grip and stable support in the shower, tub, toilet, and more. Made of corrosion-resistant stainless steel, it installs vertically, horizontally, or diagonally on tile, acrylic, fiberglass, or drywall. Includes all mounting hardware and decorative flanges for a sleek finish.

SPECIFICATIONS

What's Included:

- Vive Metal Grab Bar
- 6pc - 2" stainless steel mounting pan head screws
- 6pc - Screw Anchors

Materials:

- Stainless steel
- 2" Screws

Colors:

- Brushed Nickel, Bronze, White, Black.

Sizes:

- 12" Grab bar- 15.11" end to end when installed with decorative flanges.
- 16" Grab bar- 19.11" end to end when installed with decorative flanges.
- 24" Grab bar- 27.11" end to end when installed with decorative flanges.

Care Instructions:

- Wipe down with a damp cloth and a mild disinfectant solution or disinfectant wipe.
- Dry immediately with a lint-free cloth.

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



STABILIZING SUPPORT

Safely provides added stability in the bathroom, shower or tub.

VERSATILE POSITIONING

Secure vertically, horizontally or diagonally for optimal support in any location.



STAINLESS STEEL CONSTRUCTION

Supporting up to 440 pounds, the grab bar is made with stainless steel for exceptional durability.

STYLISH FINISH

Seamlessly blending with any decor, the metal grab bar is available in a smooth brushed nickel, white, black or bronze satin finish.



INSTALL ON MOST SURFACES

Safely install with included mounting hardware on fiberglass, acrylic, drywall and plaster walls, or on ceramic tiles with additional mounting hardware.

CONCEALED MOUNT

Includes decorative flanges to hide mounting hardware for a sleek appearance.



(239) 220-5772



sales@vivehealth.com

vive
health

SR25-7R1