



VIVE

HARD NIGHT SPLINT

SUP1035

HCPCS: L4396/L4397

Relieve foot pain and inflammation while you sleep with the Vive Hard Night Splint. Constructed with a hard composite exterior and lined with a soft, breathable padding, the night splint gently stretches the foot, maintaining a neutral position to reduce pain and inflammation related to plantar fasciitis, achilles tendonitis, and drop foot. The adjustable splint includes removable stretch wedges to elevate the foot for maximum pain relief. Use the versatile splint with or without a sock, on either the left or right foot.

SPECIFICATIONS

Sizes:

- Small: Men's: Up to 6.5 - Women's: Up to 7
- Medium: Men's: 7 - 8 - Women's: 7.5 - 9.5
- Large: Men's: 8.5 - 11 - Women's: 10 - 12.5

Materials:

- Breathable Neoprene, Strong fastener & Metal splints.

Care Instructions:

- Hand wash using cool water and mild soap.
- Dry by dabbing with towel to remove the majority of the water.
- Allow insole to air dry until all moisture is removed.

What's included:

- Single hard night splint
- 1 x spike ball
- 60 day guarantee

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



ALLEVIATE FOOT PAIN AND INFLAMMATION

Splint maintains neutral foot position while creating a gentle stretch.



CUSTOMIZE YOUR STRETCH

Use the wedges separately or combined for your desired stretch.



VERSATILE SPLINT DESIGN

Can be repositioned to fit either foot.

ALLEVIATE FOOT PAIN AND INFLAMMATION

Maintaining a neutral foot position that gently stretches the plantar fascia, the night splint reduces pain and inflammation due to plantar fasciitis, drop foot and Achilles tendonitis. Gently stretching the foot, the splint is worn throughout the night to alleviate painful first steps in the morning. The rigid exterior supports the calf, ankle and foot so you can rest comfortably.

TWO WAYS TO OPTIMIZE THE STRETCH

Customize the amount of stretch for optimal pain relief with the use of two comfortably padded stretch wedges and dual tension straps. The tension straps can be adjusted to increase angle of the foot. Made with a soft foam, the wedges elevate the foot to ensure proper positioning for maximum relief. Use each wedge separately or combined to achieve the perfect amount of stretch.

DURABLE EXTERIOR WITH A SOFT INTERIOR

Durably constructed with a hard composite material, the exterior splint provides exceptional support that can withstand extensive use. The interior of the hard night splint is comfortably padded with a breathable foam with a soft lining. The soft lining prevents skin irritation, blisters and abrasions.

VERSATILE DESIGN

Secured with three wide adjustable straps, the latex-free night splint fits both men and women. The straps are made with a strong webbing and fastener material. Each strap also has a buckle, eliminating the need to adjust the straps every night. Versatilely designed, the splint can be worn on the left or right foot.



(239) 220-5772



sales@vivehealth.com

vive
health

SR23 - BR3