



VIVE KNEE ICE WRAP

SUP1070

HCPCS: A9273

Featuring three Arctic Flex gel packs and a soft neoprene blend, the Knee Ice Wrap provides targeted hot or cold therapy with customizable compression for pain-relieving support. Each removable gel pack remains flexible when frozen to contour to the knee, reducing pain, swelling, and inflammation caused by injury, surgery, arthritis, or overuse. Secured with four adjustable straps, the versatile knee brace is reversible for use on the left or right knee. The soft knee brace is constructed with a latex-free neoprene blend that is lightweight and breathable for a comfortable fit.

SPECIFICATIONS

What's included:

- 1 Knee brace
- 1 cold packs

Materials:

- Latex-free neoprene/spandex/polyester/nylon blend
- Clay gel inside nontoxic PVC
- Extra-strength fastening material

Care Instructions:

- Remove gel packs before washing
- Machine wash in cold water with like colors
- Hang to dry

Size:

- Universal - One Size Fits Most
(Fits Right and Left Knee)

Dimensions:

- Brace: 35.25" (L) x 10" (W)
- Each gel pack: 9.5" (L) x 4.25" (W)

Colors:

- Black, Grey, Black/Blue, Grey w/Pink

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com




RECOVERY ESSENTIAL

Ice your knee with complete freedom of mobility.



COMPLETE KNEE COVERAGE

3 integrated Artic Flex ice packs cover the front, back and side of your knee.



ARTIC FLEX TECHNOLOGY

Industry leading ice packs stay cold for 30 minutes and are flexible when frozen.

SOOTHING RELIEF FOR SORE OR INJURED KNEES

Provides relief for pain and inflammation from ACL, arthritis, injuries, and surgeries.

CUSTOMIZED THERAPY

Includes three non-toxic gel packs, two in the front and one in the back for customized treatment.

HOT OR COLD RELIEF

Easy to heat or cool, enjoy up to 20 minutes of targeted temperature therapy.

ADAPTABLE TO ANY KNEE

Adjusts to fit on any adult left or right knee, up to 21" in circumference.

SECURE FIT

Stays in place without slipping for all-day knee stability and support.

WASHABLE NEOPRENE

The safe and breathable fabric will keep your skin cool and dry even with extended use.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - 10R2