



## VIVE LIFT SLING WITH OPENING

LVA2057BLU

Designed to assist in lifting, transferring, repositioning and toileting while minimizing the risk of patient or caregiver injury, the Vive Lift Sling with Opening is durably constructed with a breathable mesh material with reinforced stitching and a commode opening. The machine-washable transfer lift sling safely supports up to 400 pounds and features four reinforced sling points for use with most floor lift models and three easy-grip handles for assisting at any angle. The patient lift sling is helpful in home health care, elder and bariatric care as well as assisting the disabled and those recovering from surgery or injury.

### SPECIFICATIONS

Materials:

- Water-resistant nylon

Care and cleaning:

- Hand or machine wash on the gentle cycle in cold water with a mild detergent.
- Hang dry or machine dry on low.
- For spot cleaning, wipe the surface with a damp cloth. Allow to air dry.

Dimensions:

- 48" by 30"
- Opening: 12" by 8"
- Supports up to 400 pounds

Colors:

- Blue

What's included:

- Vive Lift Sling with Opening

**Call now!**

or send us an e-mail  
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



### CONVENIENT OPENING

Designed for easy bathroom experiences.



### ADJUSTABLE SLING HANDLES

Allows for easy repositioning for ideal angles.



### COMPATIBLE WITH MOST LIFT MODELS

Great for home health, elderly and bariatric patients.

### COMPATIBLE WITH MOST PATIENT LIFTS

Designed with four reinforced sling points to easily work with most floor lift models.

### COMMODORE OPENING

Features a 12" by 8" opening for hygienic bathroom use.

### MULTIPLE HANDLES

Includes three hand grips for easy repositioning at any angle.

### DURABLE REINFORCED DESIGN

Providing added strength and durability, the Vive Lift Sling with Opening includes reinforced corners and stitching to safely support up to 400 pounds.

### BREATHABLE MESH MATERIAL

Made with a water-resistant nylon material that is breathable for optimal patient comfort.

### MACHINE WASHABLE

Easy to clean, the Vive Lift Sling with Opening is machine washable.



(239) 220-5772



[sales@vivehealth.com](mailto:sales@vivehealth.com)

# vive health

SR24- 8R2