



VIVE

RAISED TOILET SEAT

LVA1011

HCPCS: E0244

An easy to install bathroom safety solution, the Vive raised toilet seat securely attaches to any standard or elongated toilet seat providing 5" of added height. Great for those recovery for surgery or injury, the elderly and those with limited mobility the raised seat allows individuals to safely sit or stand.

SPECIFICATIONS

What's Included:

- 1x Raised toilet seat
- 1x Lifetime Warranty

Dimensions:

- Exterior: 18" x 16.5" x 5"
- Interior: 9.5" x 8.5"
- Between the removable handles: 17"
- Fits both standard and elongated ADA compliant toilets

Materials:

- Durable composite base
- Closed-cell foam
- Corrosion-resistant aluminum

Cleaning Instructions:

- Wipe down with a damp cloth and mild soap or disinfectant.
- Air dry

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



ADDED HEIGHT TO MORE EASILY SIT OR STAND

Adding five inches of height to any standard or elongated toilet bowl, the raised toilet seat with comfortable padded handles makes it easier and safer to sit or stand without bending. A great bathroom safety solution for the elderly, those with limited mobility and those recovering from hip, knee or back surgery.



CONTOURED DESIGN WITH PADDED HANDLES

Constructed with a durable composite material that supports up to 250 pounds, the seat is contoured to provide safe and comfortable elevated seating. Removable handles are padded with a soft foam that provides a secure nonslip grip for added stability when lowering or rising from the seat.

EASY TOOL-FREE INSTALLATION ON MOST TOILETS

Installs on any standard or elongated toilet bowl without removing the existing seat. The raised seat securely attaches with an adjustable knob and rear wings for added stability.



PORTABLE, DURABLE AND EASY TO CLEAN

The durable raised toilet seat can be wiped down with regular household disinfectants or sanitizing wipes and does not interfere with regular cleaning of the toilet bowl. Easily portable, the seat riser is lightweight with removable handles.

VIVE GUARANTEE

Lifetime guarantee so you can purchase now with confidence.



(239) 220-5772



sales@vivehealth.com

vive
health

SR23-7R1