



VIVE

SUREWALK QUAD CANE

MOB1073

Introducing our SureWalk Quad Cane, a reliable mobility aid designed for comfort and stability. Crafted from sturdy yet lightweight aluminum, it supports up to 250 pounds and stands independently. The precision-made four-prong plastic base ensures enhanced stability and balance, while its height-adjustable design accommodates various heights. With an ergonomic grip handle, customizable for wrist comfort, it's the perfect companion for improved mobility and reduced strain.

SPECIFICATIONS

Materials:

➤ Aluminum.

Dimensions:

➤ 30" - 39"

➤ Base 6"-8"

Product Weight:

➤ 1.25lbs

What's included:

➤ 1x SureWalk Quad Cane

➤ 1x Allen key for adjusting grip

Colors:

➤ Black

➤ Teal

➤ Purple

Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



DURABLE ALUMINUM BUILD

Offers exceptional durability but also boasts lightweight convenience.



EASY HEIGHT ADJUSTMENTS

Height is adjustable from 30" to 39" with a simple 1" push-pin system.



ERGONOMIC HANDLE

The ergonomic grip handle effectively reduces hand and wrist fatigue.

STURDY & LIGHTWEIGHT ALUMINUM CONSTRUCTION

Crafted from heavy-duty aluminum, this quad cane offers both durability and lightweight convenience. It provides dependable support, with a weight capacity of up to 250 pounds, and stands independently when not in use.

PRECISION-MADE FOUR-PRONG PLASTIC BASE

The quad cane features a meticulously crafted plastic base that enhances stability and balance. Each cane is machine-balanced with precision, ensuring a secure footing and confident stride with every step.

HEIGHT-ADJUSTABLE DESIGN

Easily tailor the cane's height to your preference using convenient height-adjustment increments. A simple push-button mechanism allows for 1-inch height adjustments, accommodating heights between 30 inches and 39 inches as needed.

ERGONOMIC GRIP HANDLE FOR COMFORT

Experience reduced wrist strain with the adjustable grip clamp handle. Utilizing the included Allen wrench, you can customize the grip angle for an ergonomic feel, minimizing the stress on your wrist and hand. The finger indentations on the bottom of the handle ensure a natural, relaxed, and reliable grip.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - 1R2