

VIVE TENS ELECTRODES (FABRIC BACK)

- 10 SETS OF 4

RHB1033FBR40

Compatible with most TENS and EMS units, the Vive TENS unit replacement pads feature universal lead connectors. Non-irritating to the skin, this set of 52 latex and pvc-free electrode pads are pregelled for better conductivity and adhesion and are reusable, up to 25 uses per pad.

SPECIFICATIONS

Materials:

- Non-woven cloth backing
- Conductive carbon film (Low impedance)
- Conductive hydro gel
- Lead wire with 2mm DIN size or snap connector

What's included:

- 40 Replacement electrode pads with FABRIC backing, 10 sets of 4
- Lead connectors

Care Instructions:

- Device - Wipe clean with damp cloth, Wipe dry.
- Electrode pads - Wet with water and wipe clean with cloth.

Dimensions:

- 1.97" by 1.97"



Call now!

or send us an e-mail
now to find out more!



portal.vivehealth.com



(239) 220-5772



sales@vivehealth.com



ENHANCED PREGELLED PADS

Medically tested formula.



FLEXIBLE AND REUSABLE

Durably made for multiple treatment sessions.



SAFE AND EFFECTIVE

Easily replace worn out TENS electrode pads.

SET OF 40 REPLACEMENT ELECTRODE PADS

Easily replace worn out TENS unit pads with safe, FDA approved, reusable electrode pads. Conveniently packaged in sets of four pads, each set includes a resealable bag with a protective liner to provide a safe storage solution that extends the lifespan of each pad.

UNIVERSALLY COMPATIBLE CONNECTORS

Compatible with many TENS units, EMS units and other muscle stimulators, each replacement pad features a universal 2mm lead connector. The non-corrosive universal lead connectors provide current stability for safe and effective pain-relieving therapy.

SELF-ADHERING, PRE-GELLED PADS

Square shaped and backed with a non-woven cloth, the electrode pads are flexible for optimal placement on the skin. The second layer of each pad is a conductive carbon film. A premium hydrogel creates the third layer, ensuring a secure adhesion to the skin that easily releases after treatment without leaving a sticky residue behind.

LATEX-FREE REUSABLE ELECTRODES

Non-irritating to the skin, each latex and pvc-free pad has a lifespan of 20-25 uses when safely stored in the resealable pouch with the protective liner.



(239) 220-5772



sales@vivehealth.com

vive
health

SR24 - 2R1