

vive
vivehealth.com

Distributed by

vive
health

8955 Fontana Del Sol Way
Naples, FL 34109
1-800-487-3808
www.vivehealth.com

WOODEN TRANSFER BOARD

RHB1037



To see all FAQs in one place,
visit vhealth.link/f3a8f

vive
vivehealth.com

QUICK START GUIDE

CAN BE USED ON

- Wheel Chair
- Bed
- Chair



HOW TO USE

1. Place a portion of the board on the surface you are transferring to/from (1/3 suggested).

2. Support body weight on board and slowly move from one end to the other.

