

## Therapeutic 5 Zone Support Mattress

HCPCS: E0184



- High quality latex free therapeutic foam mattress provides pressure redistribution and shear/friction reduction
- Deluxe horizontal, cross-cut foam mattress provides comfort, support, and pressure redistribution over 5 therapeutic pressure zones
- Better compliance and durability than innerspring mattresses allows for improved patient satisfaction and fewer service calls
- Non skid bottom keeps mattress safely in place
- Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- Meets flammability standard CFR 1632 and CFR 1633

PRODUCTS	Description
15019	Therapeutic 5 Zone Mattress,36WX80L,1/cv
SPECIFICATIONS	Description
Warranty	Mattress: 3 Years Cover: 1 Year
Weight Capacity	350 lbs
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> ]
ACCESSORIES	Description
15019C	15019, Cover only, 1/ea, EA